

**GOAT-ANTI RABBIT IgG (H+L)
HORSE RADISH PEROXIDASE (HRP) CONJUGATED
AFFINITY PURIFIED**

| | |
|--------------------------|---|
| Catalog Number | HT-AR0710-H |
| LOT Number | 0101 |
| Target Species (Epitope) | Rabbit IgG |
| Host Species | Goat |
| Unit Size | 1 ml |
| Clonality | Polyclonal |
| Description | Whole molecule of goat anti rabbit IgG, affinity chromatography purified, HRP conjugated. |
| Specificity | The antibody reacts with the whole rabbit IgG. |
| Species Reactivity | Reacts with rabbit IgG. Reactivity with other species has not been determined. |
| Applications | Lateral Flow and ELISA. Optimal working dilutions must be determined by end user. |
| Usage | Ready to use |
| Storage/ Handling | Shipped and Store at 2-8 ° C. Upon delivery it must be aliquot in the desired volume and store. |
| Format | Liquid |
| Important Note: | During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. We recommend gently spin down the vial to dislodge any liquid in the container's cap. |

***FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC PROCEDURES.
NOT FOR HUMAN OR ANIMAL CONSUMPTION.**